Experiment 2 Verification of Newton's Second Law of Motion by Atwood Machine.

Physics Lab 1
Summer 2020-21
Department of Physics
American International University-Bangladesh

Objectives:

To establish the relationship between force and acceleration, thus verify Newton's second law of motion.

Outcomes:


After completing this experiment student should be able to answer the following questions:

- What is the relationship between force and acceleration for an object according to Newton's second law of motion?
- What is the basic concept of net force?
- How an Atwood machine can be constructed? How different forces and acceleration work for the Atwood machine.
- Why the experimental accelerations vary from the theoretical accelerations?
- What is the meaning of a linear relationship and how it looks in a graph?

Theory: Newton's Second Law of Motion


Newton's second law of motion says that FORCE equals MASS times ACCELERATION





For a particular mass:





Atwood Machine: Theory

Applying Newton's 2^{nd} Law: $\mathbf{F_{net}} = \mathbf{ma}$

• For M:

$$\mathbf{F}_{\text{net}} = \mathbf{T} - \mathbf{M}\mathbf{g} = \mathbf{M}\mathbf{a},$$

• For m:


$$\mathbf{F}_{net} = \mathbf{T} - \mathbf{mg} = \mathbf{ma}$$

• Solving these two equations: the theoretical acceleration

$$a_{th} = \frac{g}{(M+m)} (M-m)$$

• Keeping (**M**+**m**) constant at any particular place (**g** is constant), we get


$$a_{th} \propto (M-m)$$



Atwood Machine: Experiment

- M falls a distance D in time t from rest.
- Applying the knowledge of equations of motion (D=ut+1/2at²), we get the experimental acceleration

$$a_{ex} = \frac{2D}{t^2}$$



Verification of Newton's 2nd Law

- Atwood Machine: Theory Predicts

$$a_{th} \propto (M-m)$$

Atwood Machine: Experimental Result


$$\mathbf{a_{ex}} = \frac{2\mathbf{D}}{\mathbf{t}^2}$$

If we find

$$a_{ex} \propto (M - m)$$

we can say, Newton's 2nd law is verified.

 $Acceleration \propto Mass \ difference$



Apparatus


Atwood Machine:

Pulley, two hangers, different masses, string, stand and clamp.

Measurement of D and t:

Meter scale and stop watch.







Procedure:

- Hold the lighter mass on the floor attached to one end of a string. The
 heavier one attached to the other end of the string is up in the air at a
 height D from the floor. Measure D with a meter scale.
- Now release the lighter mass and measure the time the heavier mass takes to fall onto the floor. Run the experiment for 7 different mass-differences, (M m). For each run, obtain the value of the acceleration in (m/s^2) experimentally as well as theoretically. Make sure to keep total mass (M + m) always constant.
- Using Excel plot acceleration (a_{th} and a_{ex}) versus mass difference (M m) graph.

Video lecture on Procedure:



Lab Works:

Complete the data table with the calculations.

 Draw the acceleration vs mass difference graph in Excel. Plot both the accelerations (theoretical and experimental) on the same graph paper.

Analyze the result.

Discussion on Outcomes of the Lab

- What is the relationship between force and acceleration for an object according to Newton's second law of motion?
- What is the basic concept of net force?
- How an Atwood machine can be constructed? How different forces and acceleration work for the Atwood machine.
- Why the experimental accelerations vary from the theoretical accelerations?
- What is the meaning of a linear relationship and how it looks in a graph?